



ICU Patient Deterioration System

Version 1

**Clinical
User Guide**

Contents

Introduction	3
ViSIG Overview.....	3
ViSIG Index.....	3
ViSIG Alerts	3
Intended Use.....	4
Safety Information	4
Data Protection/ Privacy	5
Symbols.....	5
Trademarks	5
Licensing	6
Getting Started.....	7
Logging In.....	7
Forgotten Password	7
Users	8
User Roles	8
Dashboards	9
Patient List	10
Unit Selector	11
Patient Tiles	12
Drill Down	15
ViSIG Index Warning Thresholds	17
ViSIG Alerts	18
Silence a ViSIG Alert	19
Clear a ViSIG Silenced Alert.....	20
ViSIG Latched Alerts - Acknowledgment	20
ViSIG Alert Notice Box.....	21
Trends Charts	21
Vitals Chart.....	23
Vital Sign Warnings	24
Notifications Chart	25
Manually Input Vital Sign Data (Manual Observations)	26
Vital Sign Timeout Settings	27
Manually Admitting a Patient to a Care-Unit	28
Transferring a Patient to a different Care-Unit	29
Discharging a Patient	30
Support	31

Introduction

This guide contains all the information you need to operate ViSIG™.

ViSIG is intended for use by trained medical personnel and assumes prior knowledge of the operation of multi-parameter patient monitors.

The day-to-day administration of the system is detailed in the Administration and Configuration User Guide (Part No. 011-0814-LMAN).

ViSIG Overview

ViSIG™ is a system for assessing the risk of mortality in ICU patients. It works by examining current and past vital sign readings, as well as whether or not the patient is receiving mechanical ventilation (MV).

For each vital sign, specific patterns of readings over the previous 30 minutes that are associated with mortality are identified, summed, and along with ventilator status (room air, on MV, previously on MV), are incorporated into a score (ViSIG Index) that indicates a patient's risk of mortality.

These median measurements are then each assigned a code corresponding to the range within which the measured value lies. The four codes then represent a pattern, and the ViSIG™ system contains a set of pre-determined risk scores associated with each pattern, embodying the knowledge of which patterns of vital sign history are associated with high risk, and which are low risk.

A basic ViSIG Index is constructed by adding the scores from each vital sign. The ViSIG Index is correlated with risk of mortality, and can be used to trigger early intervention, with a visual and audible alert highlighting a deteriorating patient for clinicians to respond as protocolized.

ViSIG Index

The ViSIG Index is a single measure of a patient's current risk level for mortality.

It ranges from 0-100 and gives an up-to-date metric of a patient's probability of mortality within the ICU. As the ViSIG Index increases, so does the patient's probability of mortality without timely intervention.

The value of a patient's ViSIG Index is given a corresponding color code on the User Interface: green for baseline risk, yellow for moderately increased risk, and red for patients at highest risk of mortality.

When a ViSIG alert has activated (patient's icon is yellow or red), it means that the ViSIG Index has reached and/or surpassed the default threshold, highlighting the patient's condition to the intensivist for an appropriate response.

ViSIG Alerts

ViSIG Alerts indicate that a patient is at an increased risk of mortality. When a ViSIG alert has activated, it means that the ViSIG Index has reached and/or surpassed the default threshold, highlighting the patient's condition to the intensivist for an appropriate response.

These alerts may change over time for a patient, depending on the ViSIG Index.

Intended Use

ViSIG with alerts is intended to be used with data from already cleared sensors measuring physiological parameters, including heart rate, respiratory rate, oxygen saturation, mechanical ventilation status and blood pressure in patients being monitored in a healthcare facility.

The device provides a time series patient status index (ViSIG Index) based on patterns of change in patients monitored vital sign. The ViSIG Index is a single measure of a patient's condition and identifies those patterns in a patient's vital signs that are consistent with subsequent patient deterioration, which has been derived from pre-determined risk scores associated with each pattern, embodying the knowledge of which patterns of vital sign history are associated with high risk, and which are low risk. The ViSIG Index is based on an integrated computation evaluating changes in the parameters and their relationships to each other.

The ViSIG Index is intended for regular review by a qualified clinician. ViSIG is intended to provide additional information for use during routine patient monitoring. ViSIG is an adjunct to and is not intended to replace vital sign monitoring.

Note: Images used within this manual are provided for reference purposes only. Screens may differ based on system configuration and available parameters.

Note: Temperature is not used by the ViSIG Algorithm when calculating the ViSIG Index, however this Vital Sign can be displayed via the Vitals Chart within the Patient Tile or Drill Down Display as a configurable option.

Safety Information

Users should read and adhere to all Warnings, Cautions and Notes listed here and in the associated sections throughout this manual.

Do not use ViSIG before readings these instructions.

Warnings:

Intended Use

- ▲ ViSIG software must not be used outside of its intended use.
- ▲ ViSIG software is not for paediatric use.
- ▲ ViSIG is an adjunct to and is not intended to replace vital sign monitoring.

Alerts

- ▲ When ViSIG Alerts are silenced, there is no notification of a potentially clinically significant change in the patient's status. Observe the patient by other means when ViSIG Alerts are silenced.

Cautions:

Intended Use

- ▲ Federal (U.S.) law restricts this system to sale by or on the order of a physician.

Alerts

- ▲ Vital Sign Warnings are not intended to be used as a primary alert mechanism.
- ▲ Vital Sign Warnings are not a replacement for any primary alerting systems or procedures already in place.

General

- As a Web application, ViSIG is reliant on running via a compatible Web Browser with good network connectivity to ensure access to patient data.

Verified Web Browsers are listed in the Admin and Configuration Guide [011-0710-LMAN].

- Your system administrator will provide you access to ViSIG via a dedicated display, typically located near or at a nurse central station.
- Your system admin can also provide you access to ViSIG via a web link for your organizations computer or laptop.
- If you cannot access the system due to poor network connectivity, attempt to login via another workstation with known network connectivity.

Note: While browsers remain largely compatible between versions, full compatibility can only be guaranteed with the listed version.






Note: The various dashboards as used by ViSIG have been designed to work best on a desktop computer or laptop with a minimum screen resolution of 1024 x 768.

Data Protection/ Privacy

Clinicians and other users of the ViSIG system should be aware that, in collecting and recording patient names and data, they are responsible for complying with all applicable data protection and/or privacy laws and regulations.

Note: All patient names and data used throughout this guide are fictitious.

Symbols

Symbol	Title	Description
	Manufacturer	Indicates the medical device manufacturer
	Batch Code	Batch code so the lot or batch can be identified
	Catalogue Number	Catalogue number so product can be identified
	Consult Instructions for Use	Indicates the need for user to refer to instructions for use
	Prescription Use Only (USA only)	Device is prescription use only by a healthcare professional

Trademarks

ViSIG® is a registered name of ViSIG Health LLC.

Licensing

ViSIG uses a licensing system that provides copy protection, security, algorithm functionality and allows registered clinical users to login and view monitored patients via the various dashboards.

An activated license determines how many beds can be monitored at any given time. It is important to ensure the ViSIG system has enough licenses in place for all monitored beds.

If your system admin has not yet activated the software license, clinical users will not be able to login to view the dashboards.

If the activated license has expired, or there are more beds monitored than adequately covered by a license, then the ViSIG Index will not be calculated. A yellow border will be displayed – highlighting the ViSIG Index cannot be calculated for the additional beds.

ViSIG Demonstration Purposes Only - Hospital 1

Acuity List

Clinical User 1

Logout

CareUnit1

CareUnit2

CareUnit3

CareUnit1

BED-221022

☆Jarvis [43]
Mathew

25

1022

BED-231023

☆Gonçales [23]
Rafael

25

1023

BED-251025

☆McEdwards [M] [66]
Gavin

13

1025

BED-021002

☆Jarvis [M] [43]
Mathew

—

1002

BED-041004

Stephens [M] [62]
John

—

1004

BED-051005

☆O'Carra [M] [66]
Rónán

—

1005

BED-071007

☆Kerr [F] [71]
Melissa

—

1007

BED-211021

★Kerr [F] [43]
Laura

—

1021

CareUnit2

BED-181018

☆Lewis [F] [63]
Kim

25

1018

BED-111011

☆Edwards [M] [55]
Mark

—

1011

BED-121012

☆Bates [M] [60]
Simon

—

1012

BED-131013

☆Roberts [M] [74]
Jason

—

1013

BED-141014

Jenkins [F] [56]
Rene

—

1014

BED-151015

☆Yeoman [M] [60]
Thomas

—

1015

BED-171017

★Hockley [F] [27]
Jessica

—

1017

BED-201020

★Faherty [F] [40]
Siobhán

—

1020

CareUnit3

BED-161016

☆Bennet [45]
Sarah

13

1016

BED-191019

☆Reynolds
Debbie

13

1019

BED-091009

☆Cleeves [57]
Ann

13

1009

BED-101010

☆Sharp [62]
Lesley

13

1010

CareUnit2

BED-18

Lewis [F] [63]
Kim

25

1018

BED-11

Edwards [M] [55]
Mark

—

1011

BED-12

Bates [M] [60]
Simon

—

1012

BED-13

Roberts [M] [74]
Jason

—

1013

BED-14

Jenkins [F] [56]
René

—

1014

BED-15

Yeoman [M] [60]
Thomas

—

1015

BED-17

Hockley [F] [27]
Jessica

—

1017

BED-20

Fahey [F] [40]
Siobhán

—

1020

CareUnit3

BED-16

Bennet [45]
Sarah

13

1016

BED-19

Reynolds [61]
Debbie

13

1019

BED-09

Cleeves [57]
Ann

13

1009

BED-10

Sharp [52]
Lesley

13

1010

Note: Please inform your System Administrator if you are unable to access the system due to Licence Invalid message.

Note: Please inform your System Administrator if you notice the borders for some patient tiles displayed in yellow, indicating the existing Algorithm License requires a review as the ViSIG Index is not being computed due to a License issue.

Getting Started

Logging In

Depending on your organizations IT security policy, you will be issued login credentials.

Typically, you will receive a unique **Username** and **Password**.



Using a compatible web browser and supplied web link, navigate to the ViSIG homepage, and input your login credentials and **Login**.

- You will now be able to access the system and its various features as per the User Role assigned to you.

You can **Logout** of the system at any time by selecting the logout button, located top right of every dashboard.

Forgotten Password

If you cannot recall your Password, select **Forgot Password** link.

1. Type in your **Username or Email** and press **Submit**.
2. A password reset link will be sent to your registered email. Navigate to your email, select the reset password link, and follow the instructions as required.
3. Using your new login credentials, attempt to login again.

Note: Depending on the systems configuration the “Forgot Password” link maybe disabled. You will need to contact your System Administrator to reset your password if this is the case.

Note: A mail server will need to have been configured by your System Administrator to enable this functionality.

Users

User Roles


All system Users are allocated a **Role** which controls the level of access they have to the various functions and scope of the system as outlined Table 1: Roles / Functions Matrix.

Role	Scope
Administrator	Manage Installation, Licences, Inputs, Configurations, Users, Default Algorithm and Vitals Threshold settings. (See Part No. 011-0710-LMAN for more information)
Clinician	View and interact with the systems Clinical Dashboard and primary functionality as configured by your System Administrator.

Table 1: Roles / Functions Matrix	Roles	
	Administrator	Clinician
Functions		
System Administration:		
Manage System Licences	x	
Manage Users	x	
Manage Care Units	x	
Manage Staff assigned to Care Units	x	
Manage User Groups for Dashboard Functionality Preferences	x	
Patient Information:		
View Patient Demographic Information		x
View ViSIG Index		x
View Vital Sign Data (Vitals / Trend Data)		x
View Notifications (ViSIG Alerts / Vitals Warnings)		x
View Single or Multiple Care Units – ‘Acuity List’ or ‘My List’		x
ViSIG Alerts:		
Set Default ViSIG Alert Thresholds	x	
Silence ViSIG Alerts		x
Clear ViSIG Silenced Alerts		x
Acknowledge Latched ViSIG Alerts		x
Vital Signs:		
Manually Input Vital Sign Data		x
Set Default Vital Sign Upper / Lower Notification Thresholds	x	
Set Default Vital Sign Timeouts	x	
Apply Patient Specific Vital Signs Timeouts		x
ADT Functionality (Standalone):		
Manually Admit Patient		
Manually Transfer Patient		x
Manually Discharge Patient		x
Edit Patient Demographic Information		x

Dashboards

The Clinical Dashboard is configured according to the individual needs of the User. This section describes the ViSIG Clinical Dashboard, its main components and features and various configuration options.



The image above is of the **Clinical Dashboard** with the **Drilldown Display fixed** and the options to **display the ViSIG Trend Graph, Patient Trend Arrow, and Vital Signs** enabled within the **Patient Tiles**.

1. Unit Selector	9. Manage Patient Tab
2. Patient Information Tile	10. Drill Down Toolbar
3. Patient List (Acuity List / My List)	11. ViSIG Index, Alert Button, Patient Trend Arrow
4. ViSIG Alert Notice Box	12. ViSIG Historic Calculations Table
5. Admit Patient Button (Care Unit Specific)	13. Vital Signs Chart and Timeouts
6. Drill Down Tab	14. Logged in User, Audible Alert Settings, Patient List Options (Acuity / My List)
7. Notifications Tab	
8. Trends Tab	

Warnings	▲ ViSIG software must not be used outside of its intended use.
	▲ ViSIG software is not for paediatric use.
	▲ ViSIG is an adjunct to and is not intended to replace vital sign monitoring.

Note: As a Web application, ViSIG is reliant on running via a compatible Web Browser with good network connectivity to ensure uninterrupted access to patient data.

Note: If you cannot access the system due to poor network connectivity, contact your administrator or attempt to login via another workstation with known network connectivity.


Patient List

The **Patient List** displays all received **Patient Information** as individual **Patient Tiles** and can be configured according to the individual needs of the User. **Patient Tiles** are grouped for the various Care Units assigned to you. There are two options available for the Patient list: **Acuity List** and **My List**.

Acuity List

The **Acuity List** displays all **Patient Tiles** for a Care Unit assigned to a Clinical User.

Care Units and their associated **Patient Tiles** can be **collapsed or expanded** using the **Unit Selector** or by using the **^** or **v** for each **individual Unit** to allow you to focus on a single Care-Unit at a time, or to have an overview of all monitored patients.



Selecting a **Patient Tile** will bring up more detailed information via the **Drill Down Display**.

- **ViSIG Chart**
- **Vitals Chart**
- **Notifications Tab**
- **Trends Tab**
- **Manage Patient Tab**

Note: The Acuity List can be configured to automatically rank patients according to the **ViSIG Index**, **Admittance Time** or **Room / Bed**.

Note: A **scroll bar** will automatically appear for Users who have been assigned multiple Units with large amounts of monitored patients so that they can easily review patients and their ViSIG Index.


Note: The image shown above is of the **Clinical Dashboard** with the **Drilldown Display** set to fly in.

My List

The **My List** allows you to display all Patient Information Tiles for Patients within the **Acuity List**, you have **starred** as your favorites.

- If Patients from Multiple Care Units are assigned to a User, **they are not separated** by Care Unit as seen in the Acuity List.
- The associated Care Unit name is detailed within the Patient Tile.

You can switch between **Acuity List** and **My List** by simply selecting the appropriate list from the **Patient List** drop down, accessed by selecting your **Username**. The active list is displayed next to your Username.



Note: Patients that appear in “My List” are specific to the individual Clinical User and cannot be seen within the “My List” of another Clinical User.

Note: If enabled for your Hospital, the **Unit Selector** will always be positioned above the Patient List for both the Acuity List and My List.

Unit Selector

The “Unit Selector” displays the various **Care-Units** assigned you as a **Clinical User**.

When you first login, all assigned Care-Units within the Unit Selector will be pre-selected (**Active - White Background**), with the associated **Patient List** being expanded to display all **Patient Tiles**.

The **number of patients assigned** to the Unit, along with any **active VSI Alerts or Silenced Alerts**, and **active of silenced Vital Alerts (If enabled)** for the Unit remain visible.



Clinical Users can **activate / deactivate**, which Care-Units they want expanded or collapsed within the Patient List, by clicking on the relevant Care-Unit within the Unit Selector. (**Collapsed – Black Background**)



	Details the number of patients admitted to the Unit
	Details the number of Active ViSIG Alerts (Red) and Silenced Alert (Blue) for the Unit

Patient Tiles

The **Patient Tile** can be configured by your Administrator to display various components including:


- **Patient Information (Name, Age, Patient ID, Location).**
- The **ViSIG Index, Patient Trend Arrow, ViSIG Index Trend Graph, ViSIG Warning Notifications** for Threshold Breaches, and **Alert Status.**
- Monitored **Vital Signs** and **Warning Notifications** for Threshold Breaches.

A **white border** will be displayed for a Patient Tile selected for viewing via the **Drill Down Display.**

The **Patient Tile** will change colour according to the **ViSIG Alert Status.**

If enabled, the **Vitals Warning Icon and Alert** will only appear if a defined ViSIG or Vital Sign Threshold is breached. Depending on the ViSIG Index and alert status, the Patient Tile will appear as follows:

Example 1 – ViSIG Index (--) and associated “First ViSIG Index Calculation” Warning.



Patient Location - Bed: **ICU-4**

Hospital / Patient ID: **5241**


Patient Name: **Elaine Doe**

Age: **49**

Gender: **Female**

ViSIG Index (--) and associated “**First Calculation Icon**” indicating the system needs more vital sign data to calculate the first ViSIG Index. **Vital Sign** data displayed. **Vent Icon** displayed indicating patient receiving mechanical ventilation.

Example 2 – ViSIG Normal Threshold, Vital Sign Threshold Breach.



Patient Location - Bed: **BED-13**

Hospital / Patient ID: **1013**


Patient Name: **Jason Roberts**

Age: **74**

Gender: **Not Specified**

ViSIG Index (24), Normal Threshold Status (Green). **ViSIG Trend Graph** and **Patient Trend Arrow (Small Red)** indicating ViSIG Index has increased slightly over defined period. **Vital Sign Threshold Breach ((1) Warning Notification).**

Example 3 – ViSIG Medium Threshold Breach



Patient Location - Bed: **BED-05**

Hospital / Patient ID: **1005**


Patient Name: **Ronan O'Gara**

Age: **66**

Gender: **Male**

ViSIG Index (54), Medium Threshold Status (Yellow). **ViSIG Trend Graph** and **Vent Icon** displayed indicating patient receiving mechanical ventilation.

Example 4 – ViSIG High Threshold Breach (Active ViSIG Alert)



Patient Location - Bed: **BED-22**

Hospital / Patient ID: **1022**

Patient Name: **Mathew Jarvis**


Age: **43**

Gender: **Not Specified**

ViSIG Index (63), Alert Status (Red) and ViSIG Trend Graph. **ViSIG Alert Active for 3 Minutes.**

Patient Trend Arrow (Large Red) indicating ViSIG Index increased quickly over defined period.

Example 5 – ViSIG Silenced Alert Status



Patient Location - Bed: **BED-07**

Hospital / Patient ID: **1007**

Patient Name: **Melissa Kerr**

Age: **71**

Gender: **Female**

ViSIG Index (62), Alert / Silence Status (Blue) and VSI Trend Graph. **7 Minutes remaining before Silence Status times out.** Patient has also been **“Starred”** to appear in **“My List”**.

ViSIG Warning Notification Icons

		Red Fill – Bell Icon. Active ViSIG Alert and duration of active ViSIG Alert in HH:MM:SS .
		Red Fill - Latched Alert Icon. How long in HH:MM:SS since an Active ViSIG Alert has dropped below the High Threshold Limited, without being Silenced by a Clinical User.
		Blue Fill – Silenced Bell Icon. Silenced ViSIG Alert and the time remaining before current setting expires in XX Mins.
		ViSIG Index - - and associated First Calculation Icon. Indicating the system requires more vital sign data to calculate the first ViSIG Index. A minimum of thirty (30) minutes of data from four (4) vital signs are required to calculate the first ViSIG Index.

Vent Status and Starred Patients Icons

	Ven Icon, indicating the patient is currently receiving mechanical ventilation . Typically, the vent data is received from a configured interface engine or manually updated via the dashboard using the Drilldown Toolbar.
	Starred Patients, indicating if a patient has been selected to be viewed via the “My List” .

Vital Sign Warning Notification Icons

	Vital Warning Icon. The number of Active Vital Sign Threshold Breaches for a Patient if displayed in Patient Tile or for the Unit if displayed in the Unit Header bar.
--	--

Border Information

White Border – Patient Selected for **Drill Down Display**

Yellow Border – Algorithm License Expired or Exceeded Use Limit

Other Symbols

- -

Is displayed when there is no value to be displayed for a Vital Sign or ViSIG Index.


Sys / Dia (xx)

Rounded brackets indicate the Mean Arterial Blood Pressure (MAP) value as displayed has been received from the vitals feed. e.g., 130/77 (90)

Sys / Dia [xx]

Square brackets indicate the Mean Arterial Blood Pressure (MAP) value has been estimated using the Systolic and Diastolic values received from a vitals feed. e.g., 130/77 [90]

Example Clinical Dashboard and Patient Tiles with all states.




Note: ViSIG requires a minimum of four vital signs over a period of thirty (30) minutes to calculate the first ViSIG Index. During this time the “First Calculation” Icon displayed. Grey Marks - - will be displayed during this time or when the system is unable to calculate a ViSIG Index.

Note: Users with the authority to “Input Patient Information” will still be able to see the Patient Name and Patient ID from the **Manage Patient Tab** even if the PHI information is configured as hidden on the Patient List and Drilldown for that User.

Drill Down

The **Drill Down** display can be a **Fixed area** next to the Patient List or a **Fly In display** that is activated when a **Patient Tile** is selected. The **Drill Down display** will automatically populate when a **Patient Tile** is selected from the **Patient List**.



1. ViSIG Index and Alert Status	2. Patient Trend Arrow
3. ViSIG Index History Table	4. Drill Down Tool Bar
5. Vital Sign Chart	6. Vital Sign Timeouts









The **Drill Down** display provides additional information when reviewing a patient using the following tabs:

- **Notifications Tab**
- **Trends Tab**
- **Manage Patient Tab**

Note: ViSIG dashboards do not support zoom as part of the web application. Zooming may distort or hide certain features so display should be left at 100%.

Drill Down Tab

The Drill Down provides more information about the ViSIG Index and patient and consists of the following:

<p>○ ViSIG Index</p> <p>The ViSIG Index is a single measure of a patient's current risk level for mortality.</p> <p>It ranges from 0-100 and as the ViSIG Index increases, so does the patient's probability of mortality without timely intervention.</p>		
<p>○ Patient Trend Arrow (PTA)</p> <p>The PTA is an index, generated to measure the rate of change of the ViSIG Index over a defined period.</p> <p>A Red Arrow indicates an increase in the Index (deterioration) A Green Arrow indicates a decrease in the Index (Improvement) in the patient's condition.</p> <p>The larger the arrow, the greater the change to the Index over the defined period.</p> <p>A grey line indicates no change to the Index during the period.</p>		
		
Deterioration – Red Arrow	No Change	Improvement – Green Arrow
<p>○ ViSIG Index History Table</p> <p>Shows a list of the last five ViSIG Index calculations along with the date and time for each calculated risk index.</p>		
<p>○ Drill Down Tool Bar</p> <p>The drill down tool bar allows you carry out various functions via the main dashboard by displaying additional dialogue box for the selected function.</p>		
	Transfer Patient	Select the Transfer Patient icon to reassign the selected patient to a new location within the ViSIG System.
	Discharge Patient	Select the Discharge Patient icon to discharge the selected patient from the ViSIG System.
	Manual Observations	Select the Manual Observations icon to manually document a set of vital sign observations for a selected patient.

	Vital Sign Timeouts	Select the Vital Sign Timeout icon to manually adjust the selected patients vital sign timeouts from a configured selection to suite the patients type of monitoring. i.e. Continuous, Episodic or Hybrid.
	Vital Sign Warnings	Select the Vital Sign Warnings Icon to view the default Vital Sign Warning thresholds applied and to set patient specific thresholds if enabled.
	Close Drilldown	Select the Close Drill Down icon to return to the Patient List to review another patient as required. Only visible when the Drill Down is active.

○ **Vitals Chart**

Displays the current vital signs for the patient, their values and the date and times at which the vital signs were recorded.

- **Vital Sign Warning Notifications**

The vital sign warning notifications are displayed below the Vitals Chart

ViSIG Index Warning Thresholds

The **ViSIG Index** is **Colour Coded** according to the **Warning Threshold Ranges defined** for each monitored Care Unit. These are fixed and cannot be altered. Default set-up is outlined below:

Warning Thresholds		Threshold Warning Notification	Enable Latching	Lower Limit	Upper Limit
<input checked="" type="checkbox"/>	High Risk Range	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	59	100
<input checked="" type="checkbox"/>	Medium Risk Range			38	58
<input checked="" type="checkbox"/>	Normal Risk Range			0	37

- A **ViSIG Alert** is generated when the **lower limit of the High-Risk Range Threshold (59)** is reached or surpassed.
- ViSIG Alerts generate a **visual (Red Flashing)** and associated **audible sound** when triggered.

You can view the Care-Units **ViSIG Index Warning Thresholds** via the **Drill Down display – Manage Patient TAB**.

ViSIG Warning Thresholds			
ViSIG Warning Thresholds		Lower Limit	Upper Limit
	ViSIG High Risk Range	59	100
	ViSIG Medium Risk Range	38	58
	ViSIG Normal Risk Range	0	37




Note: Any **rapid increase in the ViSIG Index** as **highlighted by the Patient Trend Arrow**, across Normal, Medium threshold ranges, **warrants further investigation and review of the patient's status** at the bedside.

Note: How the **ViSIG Alerts** and **ViSIG Warning Thresholds** are responded to, are defined by your own Organization and Care Unit protocols.


ViSIG Alerts

All Clinical dashboards will provide a **visual alert** when the **ViSIG Index equals or exceeds** the defined **Algorithms High Risk Warning Threshold**.

- The **Unit Selector** will update to display an **Active ViSIG Alert** for a particular Care-Unit.
- The **border** of the **Patient Tile** will begin to **flash red**.
 - When an alerting Patient Tile is selected, the Drill Down Tab, is also colored red and you can then interact with the Alert Status using the **Alert Button** or **Latched Alert Button**.
 - The system can also be configured to **sound an Audible Alert** in addition to the visual alert.



The **Alert Button** and how long the **Alert Status** has been active for, will be displayed.




If **Latched Alerts** are enabled, the **Latched Alert Button** and how long the Latched Status has been active for, will be displayed.

Silence a ViSIG Alert

Silencing an ViSIG Alert acknowledges a clinical / technical response to the underlying cause of the ViSIG Alert has been carried out according to defined hospital protocol.


To **silence** an ViSIG Alert

1. Select the **Patient Information Tile** with the **Active ViSIG Alert**.
2. Using the **Drill Down** display – select the **Alert Button** from the **Drill Down Tool Bar**.



Select the **Alert Button** to Silence an active ViSIG Alert. From the drop down options available, select the time period required and **select > Silence**

3. A dialogue box will be displayed showing you a list of available time periods to Silence the ViSIG Alert for.



4. The **Patient Information Tile** border will change from an **Alert Status** to the **Silence Alert Status** indicating the ViSIG Alert has been **Silenced**.
5. The **Time Stamp** will now show time remaining, before the **Silenced Alert will time-out** e.g., 45 Minutes.

CareUnit1		CareUnit2		CareUnit3							
<div style="display: flex; justify-content: space-between;"> <div style="width: 33%;"> <p>CareUnit1</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="text-align: center;">BED-01 1001</td> <td style="text-align: center;">12</td> <td style="text-align: center;">☆</td> </tr> <tr> <td>Smith [44] Simon</td> <td></td> <td></td> </tr> </table> </div> <div style="width: 33%;"> <p style="text-align: center;">BED-02 1002</p> <div style="background-color: #00a0e3; color: white; padding: 5px; text-align: center;"> Jarvis [44] Mathew </div> <div style="display: flex; justify-content: space-around; margin-top: 5px;"> 62 19 Mins ☆ </div> </div> <div style="width: 33%;"> <p style="text-align: center;">BED-03 1003</p> <div style="background-color: #00a0e3; color: white; padding: 5px; text-align: center;"> Gonçales [24] Raúl </div> <div style="display: flex; justify-content: space-around; margin-top: 5px;"> 12 ☆ </div> </div> </div>						BED-01 1001	12	☆	Smith [44] Simon		
BED-01 1001	12	☆									
Smith [44] Simon											

Warnings


- ▲ When an ViSIG Alert is silenced, there is no additional notification of a potentially clinically significant change in the patient's status. **Observe the patient by other means when ViSIG Alerts are silenced.**

Note: The **Unit Selector** will update to show the number of Active Silenced Alerts there are for the Care Units assigned to you.

Clear a ViSIG Silenced Alert

Using the **Drill Down** display


1. Select the **ViSIG Alert Icon (Silenced)** from the **Drill Down Tool Bar**.



Select the **Silenced Button** to clear the Silenced Alert status applied.

Check the patient information and **select > Clear**.

2. A dialogue box will be displayed confirming you wish to clear the silence status.



3. Select the **Clear** button.
4. The **Silenced Alert Status will now be removed**.
5. If the ViSIG Index is still within the High Threshold Range the patient tile will start to Alert again. Apply a new Silence period as required.


Note: If the ViSIG falls **below the defined High Threshold Lower Limit value**, then the ViSIG Alert Status (both Visual and Audible) is automatically cleared from the Patient List and Drill Down Tab.

Note: If Latching is enabled, then the ViSIG Alert will remain in place until the latched alert is acknowledged.

Note: The actions of Silencing and Removing a Silenced Alert are documented via the Notifications Table.

ViSIG Latched Alerts - Acknowledgment

The system can be configured to Latch ViSIG Alerts to ensure any High Threshold Breaches of the ViSIG Index are acknowledged by a Clinical User if they subsequently drop below the threshold before the original ViSIG Alert is silenced.



Select the **Latched Alert Button** to acknowledge an ViSIG Alert status that was not silenced, before the ViSIG dropped below the High Threshold limit set. Only visible when Latched Alert status is active. Select > Confirm to remove latched alert.

Latched Alerts appear the same as an ViSIG Alert, however instead of “**Silencing**” the alert, you will need to “**Acknowledge**” the ViSIG Alert has occurred in the past.

You should always review the patients’ overall condition and risk of deterioration before acknowledging.

Once acknowledged the **Patient Tile** will revert to the current Threshold Status and colour.

ViSIG Alert Notice Box

Each Care Unit listed has an associated 'ViSIG Alert Notice Box'.



This allows caregivers to quickly identify any **ViSIG Alerts**, even if the Care Unit in question has been collapsed while reviewing patients in another Care Unit.


Trends Charts

The Trend Chart shows, in graph form, a history of the selected patient's **ViSIG** and a history of the values for each monitored **Vital Sign**.

- You can configure the system to display the various **Trend Graphs** and associated **Upper/Lower Threshold limits** for the **ViSIG** and **Monitored Vital Signs** using the **Trend Graph Configuration Icon**.
- By unselecting a vital sign or even the ViSIG from the Trend Graph Configuration settings, the graphs will be redrawn, and all unselected items will be hidden / removed.
- You can select the items again for display within the trend graph as required.
- You can temporally change the **time interval for the historic graphs** by clicking the options available.



- You can expand the trend graph to full screen using the **expand graph icon** []



Panning and zooming is available on the both the ViSIG and Vitals graphs.


You may use these features to alter the display to focus on regions of the graph. You will need a mouse and keyboard to use this functionality.

▪ **To Pan**

- Using your mouse, **left click on the graph and drag the image** (Left or Right) and release when the display shows the portion of the graph that is of interest.

▪ **To zoom**

- Using your keyboard, **hold down the Alt key** and using the mouse, **left click on the graph and drag the cursor to draw a grey box** around the area into which you would like to zoom. Repeat as required until the display shows the portion of the graph that is of interest.
- Select one of the default time intervals to reset the graph.
- Hovering your mouse pointer over individual trend graphs will reveal the values for each recorded vital for the date / time.



Note: The available time intervals displayed are defined by your Hospital Admin.

Note: Changing the time interval will cause the selected Trend Chart to be redrawn as required. The graph will also be redrawn every minute to include any additional vital sign data received.

Vitals Chart

The **Vitals Chart** can be configured to display the current vital signs for the patient, their values, and the times at which the vital signs were received.

Not all vital signs may be shown, depending on how ViSIG has been configured.

The associated vital sign data can also be color coded according to the Care Units default **Vitals Warning Thresholds**.

	Upper Limit Warning Threshold has been exceeded
	Lower Limit Warning Threshold has been exceeded
	Vital sign is between the Upper and Lower Limit Warning Thresholds

How long a Vital Sign is displayed via the Vitals Chart and used by the Algorithm when calculating the ViSIG Index, is highlighted in the Timeout box above the Vital Sign.

In this example 300 seconds or 5 Minutes.

Timeout : Standard Timeouts				
300 Secs	300 Secs	300 Secs	300 Secs	300 Secs
HR 121 17/02/2023 18:53:43	BP 130/77 (90) 17/02/2023 18:53:43	Resp -- 17/02/2023 18:48:24	Spo2 93 17/02/2023 18:53:43	Temp --

In addition, the box for each vital sign will change color to indicate whether data for the configured vital sign has ever been received in the past or stopped (timed out) for whatever reason.

	White Background – No data ever received for this vital sign.
	Grey Background – Data received in the past but has now stopped.
	Black Background – Data is currently being received for the associated vital sign.
--	No available data - Indicating that there is no valid data available, for example, if the vital signs have not been entered, or the vital sign has timed-out for consideration when calculating the current ViSIG Index.

You should always investigate the underlying causes as to why an expected vital is not being received or has stopped being transmitted to the ViSIG system. Missing vital signs can affect the sensitivity and specificity of the ViSIG Index.

More information about Timeouts and how to change them can be found in the **Vital Sign Timeout Settings** section.

Cautions	<ul style="list-style-type: none"> ▲ Vital Sign Warnings are not intended to be used as a primary alert mechanism. ▲ Vital Sign Warnings are not a replacement for any primary alerting systems or procedures already in place.
-----------------	---

Vital Sign Warnings

ViSIG can be configured to display warnings when monitored vital signs are outside of the systems defined thresholds. These warnings are indicated on the dashboard as a **visual warning icon** and **text indicator** within the **Patient List** and **Vitals Chart** as follows:

Vital Sign Warning Icon:	
Vital Sign Warning Text:	"Vital Sign" High or Low

The **values** within the Vitals Chart will be colour coded according to which threshold has been breached. The defined thresholds are displayed below the Vitals Chart for each of the configured Vital Signs.

▲	Patient Specific Value (120)	▼	Value (120)	Care-Units Default Upper Limit
▲	Patient Specific Value (40)	▼	Value (40)	Care-Units Default Lower Limit

The image below highlights Heart Rate breaching:

- The Upper Threshold (120).
- SpO2 breaching the Lower Threshold (94.0).

Timeout : Standard Timeouts				
300 Secs	300 Secs	300 Secs	300 Secs	300 Secs
HR 120 17/02/2023 19:01:35 40	BP 190/150 (163) 17/02/2023 19:01:35 91/40 (57)	Resp 25 17/02/2023 19:01:46 8	Spo2 93 17/02/2023 19:01:35 94	Temp 102.3 17/02/2023 19:02:12 95
▲ HR High			▲ Spo2 Low	

Vital Sign Warnings information will be displayed as long as the vital sign remains outside of the defined upper / lower threshold limits.

Full data as displayed is fictional

Thresholds 986532 DOE, Jane [F] [28] ICU2-BED-02

Define patient specific thresholds. Use up/down arrows to select the value required and press the save button
NOTE: This is a visual notification only. It does not replace the physiological devices or proprietary alert functionality

Heart Rate

Unit Default: 120 Patient Specific: 120
Upper: 120 Lower: 40

Blood Pressure

Blood Pressure Systolic Blood Pressure Diastolic Blood Pressure MAP

Unit Default: 190 Patient Specific: 190
Upper: 190 Lower: 91

Unit Default: 100 Patient Specific: 100
Upper: 100 Lower: 40

Unit Default: 160 Patient Specific: 160
Upper: 160 Lower: 57

Respiration Rate

Unit Default: 25 Patient Specific: 25
Upper: 25 Lower: 8

Spo2

Unit Default: N/A Patient Specific: N/A
Upper: N/A Lower: 94.0

Temperature (F)

Unit Default: 102.3 Patient Specific: 102.3
Upper: 102.3 Lower: 95.0

Save Reset Cancel

Cautions

- ▲ Vital Sign Warnings are not intended to be used as a primary alert mechanism.
- ▲ Vital Sign Warnings are not a replacement for any primary alerting systems or procedures already in place.

Notifications Chart

The Notification Chart details all historic ViSIG Alerts, Vitals Warnings and any interventions or changes made by Users.

- Using the **Drill Down** display – select the **Notifications TAB**.
 - A list of all enabled and logged Notifications will be displayed.
- Filter the Notifications Table by selecting the option required and press **OK**.
 - Algorithm Notifications
 - ViSIG Alerts, Silence Alerts and Latched Alerts
 - Acknowledged Latched Alerts
 - Cleared Silenced Status
 - ViSIG Threshold Breaches
 - Vital Notifications
 - Threshold Breach Warnings
 - Vent Notifications
 - Patient Notifications
 - Admit and Transfer

Filter Notifications

- ☒ Algorithm Notifications
- ☒ Vital Notifications
- ☒ Vent Notifications
- ☒ Patient Notifications

OK Cancel

Drill Down Notifications

ICU2-BED-02

Algorithm, Vitals, Vent, Patient

Unit	Location	Date	Type	Details	Duration	User
CareUnit1	ICU2-BED-02	17/02/2023 19:31:06	VISIG MEDIUM	50.0	00:07:59	Super Admin
CareUnit1	ICU2-BED-02	17/02/2023 19:10:34	Ventilator START		00:28:31	Super Admin
CareUnit1	ICU2-BED-02	17/02/2023 19:02:13	Vital Alert Temp LOW	36.7	00:05:11	Super Admin
CareUnit1	ICU2-BED-02	17/02/2023 19:02:00	Vital Alert Temp LOW	41.1	00:37:05	Super Admin
CareUnit1	ICU2-BED-02	17/02/2023 19:02:00	Temp LOW	106.0	00:05:15	Super Admin
CareUnit1	ICU2-BED-02	17/02/2023 19:01:36	Vital Alert Spo2 LOW	93.0	00:05:17	Super Admin
CareUnit1	ICU2-BED-02	17/02/2023 19:01:36	Vital Alert HR HIGH	121.0	00:05:17	Super Admin
CareUnit1	ICU2-BED-02	17/02/2023 19:01:36	Spo2 LOW	93.0	00:05:09	Super Admin
CareUnit1	ICU2-BED-02	17/02/2023 19:01:35	HR HIGH	121.0	00:05:10	Super Admin
CareUnit1	ICU2-BED-02	17/02/2023 19:00:49	Admit Patient			Clinical User2


< 1 >

Note: The abnormal vital warning notifications as displayed within the Notifications Tab are not intended to be used as a primary alert mechanism. It is not a replacement for any primary alerting systems or procedures already in place.


Manually Input Vital Sign Data (Manual Observations)

Vital sign data can be manually inputted for a selected patient as part of a routine observation round, or to supplement a patient monitoring device / solution that does not support that vital sign.

1. Using the **Drill Down display** – select the **Add Vitals Icon** from the **Drill Down Tool Bar**



2. A dialogue box will be displayed, with a list of Vital Signs that can be inputted.



3. Input the value for each vital sign observed, and press **Add** button.

Note: There are validation checks to ensure the value inputted does not exceed the systems expected Maximum and Minimum amounts for the inputted vital sign. You will be required to input an expected value before being able to select Add.

4. The dialogue box will disappear, and the updated Vitals Chart will now be updated to display the observed Vital Sign values and associated Timestamp.
5. The **ViSIG** will then take into consideration the manually updated vital sign when calculating the next ViSIG Index – every thirty (30) minutes.

Note: You will only be able to enter a patient's vital signs manually if your system has been configured to enable you to do so.

Note: Any manually entered vital signs will be overridden as soon as new vital sign data is received from a configured patient monitoring system that is outputting continuous vital sign data.

Vital Sign Timeout Settings

Preconfigured Vital Sign Timeouts are available to apply to patients as part of your Care-Units monitoring needs.

Timeouts help Clinical Users to identify vital signs monitored by ViSIG that have for whatever reason been disconnected or are no longer being transmitted by the primary device (Patient Monitor / Sensor) thus not being incorporated in the next calculation of the ViSIG Index.

Use the Vital Sign Timeout icon to manually adjust the selected patients vital sign timeouts from a configured selection to suite the patients type of monitoring.

Example TimeOuts:

Option	Monitoring Type	Timeouts Applied
1.	Continuous Monitoring	- All vitals set to time out after 300 Seconds (5 Minutes)
2.	Episodic Monitoring	- All vitals set to time out after 4 Hours
3.	Hybrid Monitoring (Mix of Continuous & Episodic Vital Signs)	- Some Vitals set to time out after 5 minutes - Some Vitals set to time out after 1 hour

- Using the **Drill Down display** – select the **Vital Sign Timeout Icon** from the **Drill Down Tool Bar**



- A dialogue box will be displayed, with a list of pre-configured Timeout Groups.
- Select the **appropriate Timeout option** for the monitored patient.
- Select the **Change button** to apply the new Timeout setting.
- The dialogue box will disappear and the Vitals Chart will be updated to reflect the new Timeout option chosen.

Cautions



- ▲ Vital Sign Warnings are not intended to be used as a primary alert mechanism
- ▲ Vital Sign Warnings are not a replacement for any primary alerting systems or procedures already in place

Manually Admitting a Patient to a Care-Unit

Patients should automatically be admitted to a Care-Unit as part of an configured ADT set-up.

However if you wish to **Manually Admit** a patient to a Care-Unit then action the following:

1. Select the **Admit Patient** button from the associated Care-Units Patient List.



2. The **Admit Patient** display is now shown for the **Care-Unit** selected.

3. Enter the **Patient ID, First Name, Last Name, Date of Birth (DOB), Gender, Room No# and Bed Location No#** as required and select **Save**.
4. The patient will now be added to the Care-Unit and should be visible in the **Patient List** for the Care-Unit.
5. You can now **manually add vitals** for the patient via the **Drill Down Display** to calculate the VISIG as required.

Note: If ADT is configured for the Care-Unit in question, all ADT information matching **Patient ID** and **Patient Name** received will automatically override any manually entered patient information such as **Location** and **Vital Sign Data**.

Note: Depending on the system configuration, options to Admit, Transfer and Discharge patients manually via the Dashboard may be disabled and as such not displayed via the Clinical Dashboard.

Transferring a Patient to a different Care-Unit


Patients should automatically transfer between Care-Units as part of an configured ADT set-up.

However if you wish to **Manually Transfer** a patient to a different Care-Unit then action the following:

1. Select the Patient from the **Patient List** you wish to **Transfer**.
2. Using the **Drill Down display** – select the **Transfer Patient Icon** from the **Drill Down Tool Bar**.



3. A dialogue box will be displayed, with options to **transfer a patient From and To, a new Room and Bed location**.



4. Select the Unit you wish to transfer the Patient to and input the Room and Bed location and select the **Transfer button**.
5. The patient will now be transferred to the new Care-Unit and will be visible in the **Patient List** for that Care-Unit.
 - The Notifications Table will detail the transfer request, as well as the User actioning the request.

Drill Down

Notifications

Trends


Edit Patient


ICU2

K-11


359781


DOE, John





15 mins





Algorithm, Vitals, Patient

Unit	Location	Date	Type	Details	Duration	User
ICU2	K-11	11/13/2020 14:09:22	Transfer Patient			Brett Nicolle
ICU1	IN-14	11/13/2020 11:31:17	HR HIGH	121.0	02:38:08	Super Admin
ICU1	IN-14	11/13/2020 11:31:17	SpO2 LOW	92.0	02:38:08	Super Admin
ICU1	IN-14	11/13/2020 11:31:17	Temp LOW	97.9	02:38:08	Super Admin
ICU1	IN-14	11/13/2020 11:31:17	RR HIGH	26.0	02:38:08	Super Admin
ICU1	IN-14	11/13/2020 10:54:38	Admit Patient			Brett Nicolle

Note: All historic information for the monitored patient should still be accessible providing an overview of the patient journey.


Note: Depending on the system configuration, options to Admit, Transfer and Discharge patients manually via the Dashboard may be disabled and as such not displayed via the Clinical Dashboard.

Discharging a Patient


Patients should automatically discharge from a care-unit via ADT functionality, which is configured during Installation and Set-up.

However if you wish to **manually Discharge** a patient from a Care-Unit then action the following:

1. Select the Patient from the Patient List you wish to Discharge.
2. Using the Drill Down display – select the Discharge Patient Icon from the Drill Down Tool Bar.



3. A dialogue box will be displayed, confirming you wish to discharge the selected patient.



4. Select the **Discharge Button**.
5. The patient should now be removed from the Patient List for the selected Care-Unit.
6. As the patient has now been discharged from the ViSIG system, you can no longer access the Drill Down to review the Notifications Tab and associated Graphs and Charts.

Note: If ADT receives some information after you have manually discharged the patient, the system will automatically re-admit the patient to the Unit as per the ADT Message and update as per the message received.

Note: Until the configured ADT system sends out Discharge message, the patient could re-appear in the Patient List for the Care-Unit even if you have manually Discharged the patient as ADT takes priority.

Note: Depending on the system configuration, options to Admit, Transfer and Discharge patients manually via the Dashboard may be disabled and as such not displayed via the Clinical Dashboard.

Support

Please contact your **System Administrator** for any technical issues experienced.

You can access an electronic copy of this User Guide via:

<http://www.obsmedical.com/ViSIG-user-guides/>

OBS Medical Ltd.
Unit 14, Cirencester Office Park
Tetbury Road
Cirencester
Gloucestershire
GL7 6JJ, UK
www.obsmedical.com

Doc Reference: **011-0809-LMAN-R1**
Issue Date: **24 Apr 2023**
Copyright © 2023 OBS Medical