



t+ Clinical Phone Diary Capability

Question Types


- Visual Analog Scale
- Numeric Rating Scale
- Numeric Value
- Icon Select
- Multiple Choice
- Multiple Select
- Date
- Time
- Information

Languages

- English (UK or American)
- French
- Spanish (European or American)
- Others on request


Diary Handling

- Branch logic to control sequence of questions
- Edit checks
- Automatic download of diary template
- Left to right, or right to left scales



Branch Logic

Branch logic allows the result from a multiple select question to control which question is displayed next.



Edit Checks

The diary can have unit parameters entered into it, allowing the phone to perform edit checks. If the subject enters invalid data, the phone can display a message and prompt the subject to re-enter the data within the selected range.

VAS

How much pain are you feeling

Back Next

Visual Analog Scale (VAS)

Descriptive labels at extremes of scale
 Optional tick marks
 Configurable resolution
 Optional numeric value at pointer
 Optional numbers on scale

B5-11

How much pain are you feeling?

Back Next

Numeric Rating Scale (NRS)

Scale direction can go left to right or right to left
 Configurable min and max values (max 11 boxes)

Numeric

A numeric question with upper and lower bounds and units (100.0 - 900.0)

89 Yards

Back Next

Numeric Value

Integer or decimal point control
 Optional min range
 Optional max range
 Optional "units" label

FPS

How much pain are you feeling?

Back Next

Icon Select

Up to 4 image icons of identical size can be shown
 Optional numeric ID above the icon

VRS-4

How much pain are you feeling?

No Pain
 Some pain
 Considerable pain
 Pain which could not be more severe

Back Next

Multiple Choice

Select a single answer from a list

* Branch logic can determine the next question based on the answer to this question

MultiChoice Question

A multichoice question with some initial selections

How are you feeling today?

Hello
 Hello Again
 Not today
 I'm ok

Back Uncheck Next

Multiple Select

Multiple answers can be selected

Date

Date input

Initially pre-set to current date with current date format appended

16/07/2008 DD/MM/YYYY

Back Today Next

Date

DD/MM/YYYY or MM/DD/YYYY format
 Optional format label
 Soft key to quickly enter current date (today)

Time

Time input

Initially pre-set to current time [12 hour clock] with current time format appended

10:21 AM 12 hours

Back Now Next

Time

12 or 24 hour format
 Optional format label
 Soft key to quickly enter current time (now)