



The WristCare Falls Management system comprises:

1. A WristCare Home System, consisting of everything you need for one person including a wrist unit and base unit, pre-programmed and ready for connection to the telephone line. The system is preconfigured with the following features active:
 - ✓ System alarms
 - ✓ Manual alarm
 - ✓ Non-Compliance alarm
 - ✓ Inactivity alarm
2. A Falls Prevention Night Light twin pack. These are stand alone night lights that plug into a standard 13 amp electrical socket. The nightlight can reduce falls by being a subliminal reminder to be careful.
3. The partnership and support of Vivatec in your Telecare implementation.

Explanation of alarms:

System alarms	Alarms to notify the alarm service in the case of mains power failure or disconnection, low batteries in wrist unit, low back up battery in base unit.
Manual alarm	The user has pressed the manual button on the wrist unit.
Non-Compliance alarm	The user has taken off their WristCare and cannot summon help. The alarm gives the alarm service the opportunity to prompt the user to put it back. It also provides an objective measure of compliance.
Inactivity	The user is unusually inactive and appears to be unconscious. This may be the result of a fall, or any other acute event.

Product Code: 7110

01/03/2007