



WristCare Social Services Range



Product Profile:

This new range of low cost pre-configured products meets the specific needs of Social Services.

How it supports the Preventative Technology Grant:

1. WristCare is a fully certified, award winning alarm system.
2. It extends independence and **delays admission** to Residential Care.
3. Upgradeable for joint working with Primary Care Trusts.
4. **Upgradeable** to other versions in the range.
5. WristCare supports a wide range of needs, **configurable to suit each individual**.

Alarm / Function	Manual	Falls	Wandering	Personal Safety	Assessment
Product code	7100	7110	7120	7130	7003
System	✓	✓	✓	✓	✓
Manual	✓	✓	✓	✓	✓
Non-Compliance		✓	✓	✓	✓
Inactivity		✓	✓	✓	✓
User Out			✓		
Intruder / Panic				opt	
Wellness Monitor			opt	opt	✓
Night Light		✓		opt	

System	Power failure, low batteries or other technical problem
Manual	The user has pressed the manual button on the wrist unit
Non-Compliance	The user is not wearing the wrist device and cannot summon help
Inactivity	The user is unusually inactive and appears to be unconscious
User out	The user has left their home
Intruder / Panic	The optional button by the front door has been pressed
Wellness	Objective measure of user's wellness, incl preventative alerts.
Night Light	Automatic night lights (x2), movement triggered, to reduce falls risk.
opt	Alarm or function is optional



Vivatec has been selected as a supply partner to the NHS Purchasing and Supply Agency (NHS PASA) for Telecare products and services.

📍 www.vivatec.co.uk
☎ 0870 2430 999

Vivatec
Preventative Telecare