









Patient Monitor & CO₂ Module Compatibility List

Patient Monitor	Model MIB/RS232 Interface	CO ₂ Module	CO ₂ Sensor	Verified & Supported
Philips IntelliVue Series (Serial MIB or LAN interface)	MP30	M3016A (Mainstream)	M1460A	
	MP50	M3015A/B (Microstream)	N/A	
	MP70	M3014A (Mainstream)	M2501A	
GE Carescape / Datex Series (Serial interface only)	B450	E-miniC (Airway Module)	N/A	
	B850			
	B30	N-FC-00 (Sidestream)	N/A	
	S/5	E-sCOV (Respiratory Module)	N/A	
Medtronic	CP20	N/A	N/A	

Notes:

- EA has been validated using a mix of both CO₂ Mainstream, Microstream and Sidestream capnography modules and sensors as supplied by the manufacturer of the Patient Monitors and associated accessories supplied as standard and listed here.
- Please refer to your **Patient Monitor User Guide** for more information with regards to all available accessories as required by the Patient Monitor, CO₂ Module, CO₂ Sensor and tubing. The above list is for validated & supported modules and sensors only.