

Patient Monitor & CO₂ Module Compatibility List

Patient Monitor	Model	CO ₂ Module	CO ₂ Sensor
Philips IntelliVue Series (Serial or LAN interface)	MP30	M3016A (Mainstream)	M1460A
	MP50	M3015A/B (Microstream)	N/A
	MP70	M3014A (Mainstream)	M2501A
	MX750 ³	867041 (Microstream)	N/A
GE Carescape / Datex Series (Serial interface only)	B450	E-miniC (Airway Module)	N/A
	B850		
	B30	N-FC-00 (Sidestream)	N/A
	S/5	E-sCOV (Respiratory Module)	N/A

Notes:

- EA has been validated using a mix of both CO₂ Mainstream, Microstream and Sidestream capnography modules and sensors as supplied by the manufacturer of the Patient Monitors and associated accessories supplied as standard and listed here.
- Please refer to your **Patient Monitor User Guide** for more information with regards to all available accessories as required by the Patient Monitor, CO₂ Module, CO₂ Sensor and tubing. The above list is for validated & supported modules and sensors only.
- Monitors marked with³ are compatible with both MIB and RS232 serial interface cards